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
FOR THE STANDARDIZATION OF MICROSCOPIC URINALYSIS


FAST - READ 102® system improves the standardization of microscopic urinalysis and provides precision and reproducibility such as to guarantee constant readings which will not be influenced by variations of techniques among different operators. FAST - READ 102® is made of a slide, protected by an optically transparent film, with 10 independent chambers containing a standard volume of 7 µL. After dispensing the sample on the slide application area by means of a capillary mechanism, the sediment is homogeneously distributed in the reading chamber. Each chamber is fitted with its own system for the collection of excess urine to prevent any possible contamination. Furthermore FAST - READ 102® S slide allows easy identification of the sample by using the numbers printed on the sides.

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
After completing the preparation of the urinary sediment by normal centrifugation,

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Use the pipette to mix the urinary sediment by repeatedly pressing the bulb; then

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
*HOW TO USE THE NUMBERING SYSTEM
Example: For the identification of samples*

insert a pipette into the test tube and dispense one drop of STAIN.

dispense a drop of sample onto the appropriate area of the slide.

from number 131 to 140 mark number 1 on the left side and number 3 on the right side of the slide.

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GB - FAST READ 102® CELL COUNT METHOD

Fast Read 102® is a plastic device composed of 10 counting chambers. Each chamber contains a GRID with 10 SQUARES, subdivided into 16 small squares (called SECTORS) .

After filling the counting chamber with the sample, proceed to the counting of cells distributed in N squares.

Considering that the grid contains 10 squares and each square has a dimension of 1 x 1 mm, a depth of 0.1 mm and a volume of 0.1 ul, the formula for determining the concentration of cells (cells / ml) is:

$$[\text{Cells/ml}] = \frac{(\sum \text{cells counted in N squares})}{N} \times \text{dilution factor} \times 10^4$$

10^4 = conversion from 0.1ul to 1 ml

Attention to the cells at the edges, you should only count those on two sides to avoid the risk of overestimates or underestimates.

In the example shown above we perform the cell counting on a sample concentrated 10 times:

N = 5 (number of squares considered for counting)

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The device is manufactured under license Hycor Biomedical Inc.(EU Patent Number 326349)

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CLASSIFICATION : *In order to the Directive 98/79/CE concerning medical-diagnostic in vitro devices, Biosigma SpA, defines all the components of the above-state family as medical-diagnostic in vitro devices and classifies it as "Other medical-diagnostic in vitro devices", in order to make conformable the attachment III of the directive.*

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MATERIALS: *The materials utilized for the production of the in vitro labware devices are selected on the base of the properties indicated for the individual destination use.*

The device BVS100 is manufactured in METHACRYLATE.

Methacrylate is already known as raw material for manufacturing of in vitro medical device. Its mains characteristics are: rigid, transparent, resistant to atmospheric agents, it replace the glass in every its application in which it reach high temperature (lower than 90-100 °C).

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

REF: BVS102

100 FAST-READ + 1000 PIPETTE PASTEUR in PE (BSV140) + 1 FLACONE STAIN COLORANTE

100 FAST-READ + 1000 PASTEUR PIPETTES IN PE (BSV140) + 1 BOTTLE OF STAIN DYE

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The prepared is not classified as dangerous. The Security Data Sheet is available on request.

The reagent should be stored away from direct sunlight and at a temperature between +15°C and +25°C.

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